

## Product datasheet for **RC206180L1V**

### **FIBCD1 (NM\_032843) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | FIBCD1 (NM_032843) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | FIBCD1   |
| Mammalian Cell Selection: | None   |
| Vector:                   | pLenti-C-Myc-DDK (PS100064)  |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_032843  |
| ORF Size:                 | 1383 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC206180).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_032843.3</a> , <a href="#">NP_116232.2</a>  |
| RefSeq Size:              | 3257 bp  |
| RefSeq ORF:               | 1386 bp  |
| Locus ID:                 | 84929  |
| UniProt ID:               | <a href="#">Q8N539</a>   |
| Cytogenetics:             | 9q34.12  |
| Protein Families:         | Druggable Genome, Transmembrane  |
| MW:                       | 50.7 kDa   |



[View online »](#)

**Gene Summary:**

FIBCD1 is a conserved type II transmembrane endocytic receptor that binds chitin and is located primarily in the intestinal brush border (Schlosser et al., 2009 [PubMed 19710473]). [supplied by OMIM, Apr 2010]