

Product datasheet for **RC206180L3V**

FIBCD1 (NM_032843) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FIBCD1 (NM_032843) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FIBCD1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032843
ORF Size:	1383 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206180).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032843.3 , NP_116232.2
RefSeq Size:	3257 bp
RefSeq ORF:	1386 bp
Locus ID:	84929
UniProt ID:	Q8N539
Cytogenetics:	9q34.12
Protein Families:	Druggable Genome, Transmembrane
MW:	50.7 kDa



[View online »](#)

Gene Summary:

FIBCD1 is a conserved type II transmembrane endocytic receptor that binds chitin and is located primarily in the intestinal brush border (Schlosser et al., 2009 [PubMed 19710473]). [supplied by OMIM, Apr 2010]