

Product datasheet for **RC206197L4V**


LAYN (NM_178834) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LAYN (NM_178834) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LAYN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178834
ORF Size:	1122 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206197).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178834.2
RefSeq Size:	2748 bp
RefSeq ORF:	1125 bp
Locus ID:	143903
UniProt ID:	Q6UX15
Cytogenetics:	11q23.1
Protein Families:	Transmembrane
MW:	42.3 kDa
Gene Summary:	Receptor for hyaluronate.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

[RC206197L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC206197L4V particle to overexpress human LAYN-mGFP fusion protein.