

Product datasheet for **RC206200A1V**

Human ZCCHC24 (NM_153367) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ZCCHC24 (NM_153367) AAV Particle
Tag:	Myc-DDK
Symbol:	ZCCHC24
Synonyms:	C10orf56; Z3CXXC8
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206200 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCCTGCTGTCGGCCATCGACACGAGCGCCGCCTCGGTGTACCAGCCCGCCAGCTGCTCAACTGGG
TCTACCTGTCGCTGCAGGACACGCACCAGGCTAGCGCCTTCGATGCCTTCCGGCCCGTGCCGACCGCCGG
CGCCGCACCCCGGAGCTGGCCTTCGGCAAGGGCCGCCCGAGCAGCTGGGCTCGCCCTGCACTCCAGC
TATCTCAACAGCTTCTTCCAGCTGCAGCGCGGAGAGCGCTGAGCAACAGTGTGTACAAGGGCGCCTCAC
CCTATGGCTCCCTCAACAACATCGCCGATGGCCTCAGCTCCCTCACCAGCACTTCTCAGACCTGACCCCT
CACCTCCGAGGCTCGCAAGCCAGCAAGCGGCCCCACCCAACCTACCTGTGCCACCTGTGCTTCAACAAA
GGACACTACATCAAGGACTGCCCCAGGCACGCCCCAAAGGCGAGGGCCTGACTCCATACCAGGGCAAAA
AGCGCTGCTTCGGCGAGTACAAGTGTCCCAAGTGAAGAGAAAATGGATGAGCGGGAACCTGGGCCAA
CATGGGGCAGGAGTGATCAAGTGCCACATCAACGTGTATCCACACAAGCAGAGACCCCTGGAGAAGCCC
GACGGCCTGGACGTGTCCGACCAGAGCAAGGAGCACCCGAGCACCTCTGCGAGAAGTGAAGGTCCTGG
GCTACTACTGCCGTGCGTGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206200 protein sequence Red=Cloning site Green=Tags(s) MSLLSAIDTSAASVYQPAQLLNWVYLSLQDTHQASAFDAFRPVPTAGAAPPELAFGKGRPEQLGSPLHSS YLNSFFQLQRGEALSNSVYKGASPYGSLNNIADGLSSLTEHFSDLTLTSEARKPSKRPPPNYLCHLCFNK GHYIKDCPQARPKGEGLTYPYQGKKRCFGEYKCPKCKRKWMSGNSWANMGQECIKCHINVYPHKQRPLEKP DGLDVSDQSKEHPQHLCEKCKVLGYCRRVQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_153367
ORF Size:	723 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_153367.1</u>
RefSeq Size:	4933 bp
RefSeq ORF:	726 bp
Locus ID:	219654
UniProt ID:	<u>Q8N2G6</u>
Cytogenetics:	10q22.3
MW:	26.9 kDa