

Product datasheet for **RC206228L2V**

CXorf20 (BEND2) (NM_153346) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CXorf20 (BEND2) (NM_153346) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CXorf20
Synonyms:	CXorf20
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_153346
ORF Size:	2397 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206228).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153346.2
RefSeq Size:	4720 bp
RefSeq ORF:	2400 bp
Locus ID:	139105
UniProt ID:	Q8NDZ0
Cytogenetics:	Xp22.13
MW:	87.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein which has two BEN domains in the C-terminus. These domains are found in proteins which participate in protein and DNA interactions which occur during chromatin restructuring or transcription. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011]