

Product datasheet for **RC206248L3V**

ELOVL4 (NM_022726) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ELOVL4 (NM_022726) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ELOVL4
Synonyms:	ADMD; CT118; ISQMR; SCA34; STGD2; STGD3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022726
ORF Size:	942 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206248).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022726.3
RefSeq Size:	3085 bp
RefSeq ORF:	945 bp
Locus ID:	6785
UniProt ID:	Q9GZR5
Cytogenetics:	6q14.1
Domains:	ELO
Protein Families:	Transmembrane



[View online »](#)

MW: 37.3 kDa

Gene Summary: This gene encodes a membrane-bound protein which is a member of the ELO family, proteins which participate in the biosynthesis of fatty acids. Consistent with the expression of the encoded protein in photoreceptor cells of the retina, mutations and small deletions in this gene are associated with Stargardt-like macular dystrophy (STGD3) and autosomal dominant Stargardt-like macular dystrophy (ADMD), also referred to as autosomal dominant atrophic macular degeneration. [provided by RefSeq, Jul 2008]