

Product datasheet for **RC206262L3V**

KLF3 (NM_016531) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KLF3 (NM_016531) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KLF3
Synonyms:	BKLF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016531
ORF Size:	1035 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206262).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016531.3
RefSeq Size:	2873 bp
RefSeq ORF:	1038 bp
Locus ID:	51274
UniProt ID:	P57682
Cytogenetics:	4p14
Domains:	zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

MW: 38.6 kDa

Gene Summary: Binds to the CACCC box of erythroid cell-expressed genes. May play a role in hematopoiesis (By similarity).[UniProtKB/Swiss-Prot Function]