

## Product datasheet for RC206322L3

### Fuc-TIX (FUT9) (NM\_006581) Human Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Fuc-TIX (FUT9) (NM_006581) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Fuc-TIX
Synonyms:	Fuc-TIX
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206322).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_006581
ORF Size:	1077 bp



[View online >](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_006581.2</a>
<b>RefSeq Size:</b>	12815 bp
<b>RefSeq ORF:</b>	1080 bp
<b>Locus ID:</b>	10690
<b>UniProt ID:</b>	<a href="#">Q9Y231</a>
<b>Cytogenetics:</b>	6q16.1
<b>Domains:</b>	Glyco_transf_10
<b>Protein Families:</b>	Transmembrane
<b>Protein Pathways:</b>	Glycosphingolipid biosynthesis - globo series, Glycosphingolipid biosynthesis - lacto and neolacto series, Metabolic pathways
<b>MW:</b>	42 kDa
<b>Gene Summary:</b>	The protein encoded by this gene belongs to the glycosyltransferase family. It is localized to the golgi, and catalyzes the last step in the biosynthesis of Lewis X (LeX) antigen, the addition of a fucose to precursor polysaccharides. This protein is one of the few fucosyltransferases that synthesizes the LeX oligosaccharide (CD15) expressed in the organ buds progressing in mesenchyma during embryogenesis. It is also responsible for the expression of CD15 in mature granulocytes. A common haplotype of this gene has also been associated with susceptibility to placental malaria infection. [provided by RefSeq, Nov 2011]

Product images:



Circular map for RC206322L3