

## Product datasheet for **RC206332L4V**

### ACBD7 (NM\_001039844) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ACBD7 (NM_001039844) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ACBD7
Synonyms:	bA455B2.2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001039844
ORF Size:	264 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206332).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001039844.2</a>
RefSeq Size:	3383 bp
RefSeq ORF:	267 bp
Locus ID:	414149
UniProt ID:	<a href="#">Q8N6N7</a>
Cytogenetics:	10p13
MW:	9.8 kDa
Gene Summary:	Binds medium- and long-chain acyl-CoA esters.[UniProtKB/Swiss-Prot Function]


[View online »](#)