

## Product datasheet for **RC206339L4V**

### **TUBB4B (NM\_006088) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | TUBB4B (NM_006088) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | TUBB4B   |
| Synonyms:                 | Beta2; LCAEOD; TUBB2; TUBB2C   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_006088  |
| ORF Size:                 | 1335 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC206339).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_006088.5</a>  |
| RefSeq Size:              | 1591 bp  |
| RefSeq ORF:               | 1338 bp  |
| Locus ID:                 | 10383  |
| UniProt ID:               | <a href="#">P68371</a>   |
| Cytogenetics:             | 9q34.3   |
| Domains:                  | tubulin  |
| Protein Families:         | Druggable Genome   |



[View online »](#)

**Protein Pathways:** Gap junction, Pathogenic Escherichia coli infection

**MW:** 49.8 kDa

**Gene Summary:** Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain. [UniProtKB/Swiss-Prot Function]