

Product datasheet for **RC206378L4V**

BPIL1 (BPIFB2) (NM_025227) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BPIL1 (BPIFB2) (NM_025227) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BPIL1
Synonyms:	BPIL1; C20orf184; dj726C3.2; LPLUNC2; RYSR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025227
ORF Size:	1374 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206378).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025227.1
RefSeq Size:	1936 bp
RefSeq ORF:	1377 bp
Locus ID:	80341
UniProt ID:	Q8N4F0
Cytogenetics:	20q11.21
Domains:	BPI1
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane


[View online »](#)

MW: 49.1 kDa

Gene Summary: This gene encodes a member of the lipid transfer/lipopolysaccharide binding protein (LT/LBP) gene family. It is highly expressed in hypertrophic tonsils. This gene and three other members of the LT/LBP gene family form a cluster on the long arm of chromosome 20.
[provided by RefSeq, Jul 2008]