

## Product datasheet for **RC206392L4V**

### FGFBP2 (NM\_031950) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	FGFBP2 (NM_031950) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FGFBP2
Synonyms:	HBP17RP; KSP37
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031950
ORF Size:	669 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206392).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031950.2</a> , <a href="#">NP_114156.1</a>
RefSeq Size:	1188 bp
RefSeq ORF:	672 bp
Locus ID:	83888
UniProt ID:	<a href="#">Q9BYJ0</a>
Cytogenetics:	4p15.32
Protein Families:	Secreted Protein
MW:	24.6 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the fibroblast growth factor binding protein family. The encoded protein is a serum protein that is selectively secreted by cytotoxic lymphocytes and may be involved in cytotoxic lymphocyte-mediated immunity. An increase in the amount of gene product may be associated with atopic asthma and mild extrinsic asthma.[provided by RefSeq Staff, Oct 2008]