

Product datasheet for **RC206402A1V**

Human HDHD3 (NM_031219) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HDHD3 (NM_031219) AAV Particle
Tag:	Myc-DDK
Symbol:	HDHD3
Synonyms:	2810435D12Rik; C9orf158
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206402 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCACACCGGCTGCAGATACGACTGCTGACGTGGGATGTGAAGGACACGCTGCTCAGGCTCCGCCACC
CCTTAGGGGAGGCCTATGCCACCAAGGCCCGGGCCATGGGCTGGAGGTGGAGCCCTCAGCCCTGGAACA
AGGCTTCAGGCAGGCATACAGGGCTCAGAGCCACAGCTTCCCCAACTACGGCCTGAGCCACGGCCTAACC
TCCCGCCAGTGGTGGCTGGATGTGGTCCTGCAGACCTCCACCTGGCGGGTGTCCAGGATGCTCAGGCTG
TAGCCCCATCGCTGAACAGCTTTATAAAGACTTCAGCCACCCCTGCACCTGGCAGGTGTTGGATGGGGC
TGAGGACACCCTGAGGGAGTGCCGCACACGGGGTCTGAGACTGGCAGTGATCTCCAACCTTGACCGACGG
CTAGAGGGCATCCTGGGGGGCCTTGGCCTGCGTGAACACTTCGACTTTGTGCTGACCTCCGAGGCTGCTG
GCTGGCCCAAGCCGACCCCGCATTTTCCAGGAGGCCTTGGGCTTGCTCATATGGAACAGTAGTGGC
AGCCCATGTTGGGGATAATTACCTCTGCGATTACCAGGGGCCTCGGGCTGTGGGCATGCACAGCTTCCTG
GTGGTTGGCCACAGGCACTGGACCCCGTGGTCAGGGATTCTGTACCTAAAGAACACATCTCCCTCTC
TGGCCCATCTCTGCCTGCCCTGACTGCCTAGAGGGCTCAACTCCAGGGCTT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206402 protein sequence Red=Cloning site Green=Tags(s) MAHRLQIRLLTWDVKDTLLRLRHPLGEAYATKARAHGLEVEPSALEQGFRQAYRAQSHSFPNYGLSHGLT SRQWWLDVVLQTFHLAGVQDAQAVAPIAEQLYKDFSHPCTWQVLDGAEDTLRECRTRGLRLAVISNFDRR LEGILGGLGLREHFDVLTSEAAGWPKPDPRIFQEALRLAHMEPVVAHVGDNYLCDYQGPRAVGMHSFL VVGPAQALDPVVRDSVPKEHILPSLAHLLPALDCLEGSTPGL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_031219
ORF Size:	753 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_031219.2</u>
RefSeq Size:	1787 bp
RefSeq ORF:	756 bp
Locus ID:	81932
UniProt ID:	<u>Q9BSH5</u>
Cytogenetics:	9q32
MW:	28 kDa