


Product datasheet for **RC206402L4V**

HDHD3 (NM_031219) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HDHD3 (NM_031219) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HDHD3
Synonyms:	2810435D12Rik; C9orf158
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031219
ORF Size:	753 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206402).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031219.2
RefSeq Size:	1787 bp
RefSeq ORF:	756 bp
Locus ID:	81932
UniProt ID:	Q9BSH5
Cytogenetics:	9q32
Domains:	Hydrolase
MW:	28 kDa


[View online »](#)

Product images:

[RC206402L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC206402L4V particle to overexpress human HDHD3-mGFP fusion protein.