

## Product datasheet for **RC206404L3V**

### **Bcl G (BCL2L14) (NM\_138722) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bcl G (BCL2L14) (NM_138722) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Bcl G
Synonyms:	BCLG
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_138722
ORF Size:	981 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206404).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138722.1</a> , <a href="#">NP_620048.1</a>
RefSeq Size:	1886 bp
RefSeq ORF:	984 bp
Locus ID:	79370
UniProt ID:	<a href="#">Q9BZR8</a>
Cytogenetics:	12p13.2
Protein Families:	Druggable Genome
MW:	36.6 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. Overexpression of this gene has been shown to induce apoptosis in cells. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported for this gene. [provided by RefSeq, May 2009]