

Product datasheet for **RC206431L4V**

AXL (NM_021913) Human Tagged ORF Clone Lentiviral Particle


Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | AXL (NM_021913) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | AXL |
| Synonyms: | ARK; JTK11; Tyro7; UFO |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_021913 |
| ORF Size: | 2682 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC206431). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_021913.2 |
| RefSeq Size: | 5014 bp |
| RefSeq ORF: | 2685 bp |
| Locus ID: | 558 |
| UniProt ID: | P30530 |
| Cytogenetics: | 19q13.2 |
| Protein Families: | Druggable Genome, Protein Kinase, Transmembrane |
| MW: | 98.2 kDa |


[View online »](#)

Gene Summary:

The protein encoded by this gene is a member of the Tyro3-Axl-Mer (TAM) receptor tyrosine kinase subfamily. The encoded protein possesses an extracellular domain which is composed of two immunoglobulin-like motifs at the N-terminal, followed by two fibronectin type-III motifs. It transduces signals from the extracellular matrix into the cytoplasm by binding to the vitamin K-dependent protein growth arrest-specific 6 (Gas6). This gene may be involved in several cellular functions including growth, migration, aggregation and anti-inflammation in multiple cell types. Alternative splicing results in multiple transcript variants of this gene. [provided by RefSeq, Jul 2013]

Product images:


[RC206431L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC206431L4V particle to overexpress human AXL-mGFP fusion protein.