

## Product datasheet for **RC206435L3V**

### Alpha 2 Antiplasmin (SERPINF2) (NM\_000934) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Alpha 2 Antiplasmin (SERPINF2) (NM_000934) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Alpha 2 Antiplasmin
Synonyms:	A2AP; AAP; ALPHA-2-PI; alpha2AP; API; PI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000934
ORF Size:	1473 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206435).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000934.3</a>
RefSeq Size:	2284 bp
RefSeq ORF:	1476 bp
Locus ID:	5345
UniProt ID:	<a href="#">P08697</a>
Cytogenetics:	17p13.3
Domains:	SERPIN
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

**Protein Pathways:** Complement and coagulation cascades

**MW:** 54.6 kDa

**Gene Summary:** This gene encodes a member of the serpin family of serine protease inhibitors. The protein is a major inhibitor of plasmin, which degrades fibrin and various other proteins. Consequently, the proper function of this gene has a major role in regulating the blood clotting pathway. Mutations in this gene result in alpha-2-plasmin inhibitor deficiency, which is characterized by severe hemorrhagic diathesis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]