

Product datasheet for **RC206458A1V**

Human SPINLW1 (EPPIN) (NM_020398) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SPINLW1 (EPPIN) (NM_020398) AAV Particle
Tag: Myc-DDK
Symbol: SPINLW1
Synonyms: CT71; CT72; dJ461P17.2; SPINLW1; WAP7; WFDC7
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC206458 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGATCTTCTGGACTTTTGAGCCTCCTGGTCTATTCGTCCTCTAGCGAATGTCCAGGGACCTGGTC
 TGACTGATTGGTTATTTCCAGGAGATGTCCAAAATCAGAGAAGAATGTGAATCCAAGAAAGGATGT
 GTGTACAAAGGACAGACAATGCCAGGACAACAAGAAGTGTGTCTTCAGCTGCGGAAAAAATGTTTA
 GATCTCAACAAGATGTATGCGAAATGCCAAAAGAACTGGCCCTGCCTGGCTTATTTCTTCATTGGT
 GGTATGACAAGAAAGATAATACTTGTCCATGTTTGTCTATGGTGGCTGCCAGGGAACAATAACAACCT
 CCAATCCAAGCCAACCTGCCTGAACACCTGCAAGAATAAACGCTTTCCC

A**CGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206458 protein sequence
 Red=Cloning site Green=Tags(s)

MGSSGLLSLLVLFVLLANVQGPGLTDWLFPRRCPKIREECFQERDVCTKDRQCQDNKKCCVFSCGKKCL
 DLKQDVCEMPKETGPCLAYFLHWWYDKDNTCSMFVYGGCQGNNNFQSKANCLNTCKNKRFP

TR**TRPLEQKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_020398


[View online »](#)

ORF Size:	399 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_020398.2</u>
RefSeq Size:	1998 bp
RefSeq ORF:	402 bp
Locus ID:	57119
UniProt ID:	<u>O95925</u>
Cytogenetics:	20q13.12
MW:	15.3 kDa