

Product datasheet for **RC206462A1V**

Human SAP30BP (NM_013260) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SAP30BP (NM_013260) AAV Particle
Tag:	Myc-DDK
Symbol:	SAP30BP
Synonyms:	HCNGP; HTRG; HTRP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206462 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGGGAAGAAGAATGTTCTGTCGTCCTCTCGCAGTTTACGCGGAAGATTCAGAGCCCGAGTCTGATG
GCGAGGCTGGAATCGAGGCGGTGGGCAGCGCGGCTGAGGAGAAAGCGGATTGGTATCTGATGCCTATGG
GGAGGATGACTTTTCTCGTCTAGGGGTGATGAAGATGGTTATGAAGAAGAAGAAGATGAGAACAGTAGA
CAGTCGGAAGATGACGATTCAGAGACTGAAAAACCTGAGGCTGATGACCCAAAGGATAATACAGAAGCAG
AAAAGCGAGACCCAGGAACCTGTGGCTCCTTTTCTGAAAGAGTTCGGAACATGTCGCCTGATGAAAT
CAAGATCCCGCCAGAACCCCTGGCAGATGTTCAAATCACTTGCAAGACAAGATCCAGAAGCTTTATGAA
CGAAAGATAAAGGAGGAATGGATATGAACTACATTATCCAAAGGAAGAAAGAAATTTGGAACCCCTAGCA
TCTACGAGAAGCTGATCCAGTTCTGTGCCATTGACGAGCTTGGCACCAACTACCCAAAGGATATGTTTGA
TCCCATGGCTGGTCTGAGGACTCCTACTATGAGGCATTAGCCAAGGCCAGAAAATTGAGATGGACAAA
TTGAAAAGGCCAAAAAGGAGCGAACAAAAATTGAGTTTGTGACGGGCACCAAAAAGGCCACACGACCA
ACGCCACGTCCACCACCACTACCACTGCCAGCACAGCTGTTGCAGATGCTCAGAAGAGAAAGAGCAAGTG
GGATTCGGCTATCCAGTGACAACGATAGCCCAGCCCACCATCCTCACCACCACAGCCACCCTGCCAGCT
GTTGTACGGTCACCACCAGCGCCAGCGCTCCAAGACCACCGTCATCTCTGCTGTGGGCACCATTTGTA
AGAAGGCCAAGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206462 protein sequence Red=Cloning site Green=Tags(s) MAGKKNVLSLAVYAEDSEPESDGEAGIEAVGSAAEEKGGLVSDAYGEDDFSRLLGGDEGGYEEEDENSR QSEDDSETKEPEADDPKDNTEAEKRDPELVASFSEVRNMSPDEIKIPPEPPGRCSNHLQDKIQKLYE RKIKEGMDMNYIIQRKKEFRNPSIYEKLIQFCAIDELGTNYPKDMFDPHGWSSEDSYYEALAKAQKIEMDK LEKAKKERTKIEFVTGKKGTTTATSTTTTASTAVADAQKRKSKWDSAIPVTITIAQPTILTTTATLPA VVTVTTSASGSKTTVISAVGTIVKKAKQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_013260
ORF Size:	924 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_013260.5</u>
RefSeq Size:	2758 bp
RefSeq ORF:	927 bp
Locus ID:	29115
UniProt ID:	<u>Q9UHR5</u>
Cytogenetics:	17q25.1
MW:	33.9 kDa