

Product datasheet for **RC206462L4V**

SAP30BP (NM_013260) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SAP30BP (NM_013260) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SAP30BP
Synonyms:	HCNGP; HTRG; HTRP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013260
ORF Size:	924 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206462).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013260.5
RefSeq Size:	2758 bp
RefSeq ORF:	927 bp
Locus ID:	29115
UniProt ID:	Q9UHR5
Cytogenetics:	17q25.1
MW:	33.9 kDa



[View online »](#)

Gene Summary:

Induces cell death. May act as a transcriptional corepressor of a gene related to cell survival. May be involved in the regulation of beta-2-microglobulin genes.[UniProtKB/Swiss-Prot Function]