

Product datasheet for RC206495L1

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Granzyme B (GZMB) (NM_004131) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Granzyme B (GZMB) (NM_004131) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: Granzyme B

Synonyms: C11; CCPI; CGL-1; CGL1; CSP-B; CSPB; CTLA1; CTSGL1; HLP; SECT

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide


The ORF insert of this clone is exactly the same as(RC206495).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_004131

ORF Size: 741 bp



Granzyme B (GZMB) (NM_004131) Human Tagged Lenti ORF Clone - RC206495L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 004131.3</u>

 RefSeq Size:
 941 bp

 RefSeq ORF:
 744 bp

 Locus ID:
 3002

 UniProt ID:
 P10144

 Cytogenetics:
 14q12

Domains: Tryp SPc

Protein Families: Druggable Genome, Protease

Protein Pathways: Allograft rejection, Autoimmune thyroid disease, Graft-versus-host disease, Natural killer cell

mediated cytotoxicity, Type I diabetes mellitus

MW: 27.7 kDa

Gene Summary: This gene encodes a member of the granzyme subfamily of proteins, part of the peptidase S1

family of serine proteases. The encoded preproprotein is secreted by natural killer (NK) cells and cytotoxic T lymphocytes (CTLs) and proteolytically processed to generate the active protease, which induces target cell apoptosis. This protein also processes cytokines and


degrades extracellular matrix proteins, and these roles are implicated in chronic

inflammation and wound healing. Expression of this gene may be elevated in human patients


with cardiac fibrosis. [provided by RefSeq, Sep 2016]



Product images:



Circular map for RC206495L1



Double digestion of RC206495L1 using Sgfl and Mlul