

Product datasheet for **RC206496A1V**

Human CISH (NM_145071) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CISH (NM_145071) AAV Particle
Tag:	Myc-DDK
Symbol:	CISH
Synonyms:	BACTS2; CIS; CIS-1; G18; SOCS
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206496 representing NM_145071 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTCCTCTGCGTTCAGGGACCTCGTCCTTTGCTGGCTGTGGAGCGGACTGGGCAGCGGCCCTGTGGG
CCCCGTCCTGGAAGTCCCAAGCCAGTCATGCAGCCCTTGCCTGCTGGGGCCTTCCTCGAGGAGGTGGC
AGAGGGTACCCAGCCAGACAGAGAGTGAGCCAAAGGTGCTGGACCCAGAGGAGGATCTGCTGTGCATA
GCCAAGACCTTCTCTACCTTCGGGAATCTGGCTGGTATTGGGGTTCCATTACGGCCAGCGAGGCCCGAC
AACACCTGCAGAAGATGCCAGAAGGCACGTTCTTAGTACGTGACAGCACGCCCCAGCTACCTGTTTAC
GCTGTCACTGAAAACCACTCGTGGCCCCACCAATGTACGCATTGAGTATGCTGACTCCAGCTTCCGTCTG
GACTCCAAGTCTGTCCAGGCCACGCATCCTGGCCTTTCCGGATGTGGTCAGCCTTGTGCAGCACTATG
TGGCCTCCTGCACTGCTGATACCCGAAGCGACAGCCCCGATCCTGCTCCACCCCGGCCCTGCCTATGCC
TAAGGAGGATGCGCCTAGTGACCCAGCACTGCCTGCTCCTCCACAGCCACTGCTGTACACCTAAACTG
GTGCAGCCCTTTGTACGAGAAGCAGCGCCCGCAGCCTGCAACACCTGTGCCGCTTGTATCAACCGTC
TGGTGGCCGACGTGGACTGCCTGCCACTGCCCCGGCGCATGGCCGACTACCTCCGACAGTACCCCTTCCA
GCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206496 representing NM_145071 Red=Cloning site Green=Tags(s) MVLCVQGPRLLAVERTGQRPLWAPSLELPKPVMQPLPAGAFLEEVAEGTPAQTESEPKVLDPEEDLLCI AKTF SYLRESGWYGSITASEARQHLQKMPEGTFLVRDSTHPSYLF TLSVKTTRGPTNVRIEYADSSFRL DSNCLSRPRILAFPDVSVLVQHYVASCTADTRSDSPDPAPTPALPMPKEDAPSDPALPAPPATAVHLKL VQPFVRRSSARSLQHLCLVINRLVADVDCPLPRRMADYLRQYPFQL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_145071
ORF Size:	774 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_145071.1</u>
RefSeq Size:	2035 bp
RefSeq ORF:	777 bp
Locus ID:	1154
UniProt ID:	<u>Q9NSE2</u>
Cytogenetics:	3p21.2
MW:	28.5 kDa