

Product datasheet for **RC206497L1V**

Granzyme M (GZMM) (NM_005317) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Granzyme M (GZMM) (NM_005317) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Granzyme M
Synonyms:	LMET1; MET1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005317
ORF Size:	771 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206497).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005317.2
RefSeq Size:	958 bp
RefSeq ORF:	774 bp
Locus ID:	3004
UniProt ID:	P51124
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane
MW:	27.4 kDa



[View online »](#)

Gene Summary:

Human natural killer (NK) cells and activated lymphocytes express and store a distinct subset of neutral serine proteases together with proteoglycans and other immune effector molecules in large cytoplasmic granules. These serine proteases are collectively termed granzymes and include 4 distinct gene products: granzyme A, granzyme B, granzyme H, and the protein encoded by this gene, granzyme M. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]