

Product datasheet for RC206503L1

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Carbonic Anhydrase XIV (CA14) (NM_012113) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Carbonic Anhydrase XIV (CA14) (NM_012113) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: Carbonic Anhydrase XIV

Synonyms: CAXiV Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC206503).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_012113

ORF Size: 1011 bp





Carbonic Anhydrase XIV (CA14) (NM_012113) Human Tagged Lenti ORF Clone - RC206503L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 012113.1</u>

 RefSeq Size:
 1757 bp

 RefSeq ORF:
 1014 bp

 Locus ID:
 23632

 UniProt ID:
 Q9ULX7

Cytogenetics: 1q21.2

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Nitrogen metabolism

MW: 37.7 kDa


Gene Summary: Carbonic anhydrases (CAs) are a large family of zinc metalloenzymes that catalyze the

reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. They show extensive diversity in tissue distribution and in their subcellular localization. CA XIV is predicted to be a type I membrane protein and shares highest sequence similarity with the other transmembrane CA isoform, CA XII; however, they have different patterns of tissue-specific expression and thus


may play different physiologic roles. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC206503L1



Double digestion of RC206503L1 using Sgfl and Mlul