

Product datasheet for **RC206520A1V**

Human CGREF1 (NM_006569) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CGREF1 (NM_006569) AAV Particle
Tag:	Myc-DDK
Symbol:	CGREF1
Synonyms:	CGR11
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206520 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTTACCTTTGACGATGACAGTGTTAATCCTGCTGCTGCTCCCCACGGGTGAGGCTGCCCCAAAGGATG
GAGTCACAAGGCCAGACTCTGAAGTGCAGCATCAGCTCCTGCCAACCCCTTCCAGCCAGGCCAGGAGCA
GCTCGGACTTCTGCAGAGCTACCTAAAGGGACTAGGAAGGACAGAAGTGCAACTGGAGCATCTGAGCCGG
GAGCAGGTTCTCCTCTACCTCTTTGCCCTCCATGACTATGACCAGAGTGGACAGCTGGATGGCCTGGAGC
TGCTGTCCATGTTGACAGCTGCTTGGCCCTGGAGCTGCCAACTCTCCTACCACCAACCCGGTGATATT
GATAGTGGACAAAGTGCTCGAGACGCAGGACCTGAATGGGGATGGGCTCATGACCCTGCTGAGCTCATC
AACTTCCCGGGAGTAGCCCTCAGGCACGTGGAGCCCGGAGAGCCCTTGCTCCATCTCCTCAGGAGCCAC
AAGCTGTTGGAAGGCAGTCCCTATTAGCTAAAAGCCCATTAAAGACAAGAAACACAGGAAGCCCTGGTCC
CAGAGAAGAAGCAAAGGCCAGGTAGAGGCCAGAAGGGAGTCTTTGGATCCTGTCCAGGAGCCTGGGGGC
CAGGCAGAGGCTGATGGAGATGTTCCAGGGCCAGAGGGGAAGCTGAGGGCCAGGCAGAGGCTAAAGGAG
ATGCCCTGGGGCCAGAGGGGAAGCTGGGGCCAGGCAGAGGCTGAAGGAGATGCCCCGGGGCCAGAGG
GGAAGCTGGGGCCAGGCAGAGGCCAGGAGAATGGAGAGGAGCCAAGGAATTCAGGGGAAACACTG
GAGTCTAAGAACACCCAAAATGACTTTGAGGTGCACATTGTTCAAGTGGAGAATGATGAGATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206520 protein sequence Red=Cloning site Green=Tags(s) MLPLTMTVLILLLLPTGQAAPKDGVTDPSEVQHQLLPNPFQPGQEQLGLLQSYLKGLGRTEVQLEHLSEQVLLYLFALHDYDQSGQLDGLELLSMLTAALAPGAANSPTTNPVILIVDKVLETQDLNGDGLMTPAELINFPGVALRHVEPGEPLAPSPQEPQAVGRQSLLAKSPLRQETQEAPGPREEAKGQVEARRESLDPVQEPGGQAEADGDVPGPRGEAEQGAEAKGDAPGPRGEAGGQAEAGDAPGPRGEAGGQAEARENGEEAKELPGETLESKNTQNDFEVHIVQVENDEI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006569
ORF Size:	903 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006569.2
RefSeq Size:	1934 bp
RefSeq ORF:	957 bp
Locus ID:	10669
UniProt ID:	Q99674
Cytogenetics:	2p23.3
MW:	31.9 kDa