

## Product datasheet for **RC206536L4V**

### **P2X6 (P2RX6) (NM\_005446) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | P2X6 (P2RX6) (NM_005446) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | P2X6   |
| Synonyms:                 | P2RXL1; P2X6; P2XM   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_005446  |
| ORF Size:                 | 1323 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC206536).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_005446.3</a>  |
| RefSeq Size:              | 2754 bp  |
| RefSeq ORF:               | 1326 bp  |
| Locus ID:                 | 9127   |
| UniProt ID:               | <a href="#">O15547</a>   |
| Cytogenetics:             | 22q11.21   |
| Protein Families:         | Druggable Genome, Ion Channels: ATP Receptors, Transmembrane   |
| Protein Pathways:         | Calcium signaling pathway, Neuroactive ligand-receptor interaction   |



[View online »](#)

**MW:** 48.6 kDa

**Gene Summary:** The protein encoded by this gene belongs to the family of P2X receptors, which are ATP-gated ion channels and mediate rapid and selective permeability to cations. This gene is predominantly expressed in skeletal muscle, and regulated by p53. The encoded protein is associated with VE-cadherin at the adherens junctions of human umbilical vein endothelial cells. Alternative splicing results in multiple transcript variants. A related pseudogene, which is also located on chromosome 22, has been identified. [provided by RefSeq, Apr 2009]