

## Product datasheet for **RC206541L3V**

### ROR gamma (RORC) (NM\_005060) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ROR gamma (RORC) (NM_005060) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ROR gamma
Synonyms:	IMD42; NR1F3; RORG; RZR-GAMMA; RZRG; TOR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005060
ORF Size:	1554 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206541).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005060.3</a>
RefSeq Size:	3084 bp
RefSeq ORF:	1557 bp
Locus ID:	6097
UniProt ID:	<a href="#">P51449</a>
Cytogenetics:	1q21.3
Protein Families:	Druggable Genome, Transcription Factors
MW:	58.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a DNA-binding transcription factor and is a member of the NR1 subfamily of nuclear hormone receptors. The specific functions of this protein are not known; however, studies of a similar gene in mice have shown that this gene may be essential for lymphoid organogenesis and may play an important regulatory role in thymopoiesis. In addition, studies in mice suggest that the protein encoded by this gene may inhibit the expression of Fas ligand and IL2. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]