

Product datasheet for **RC206543L3V**

STOML1 (NM_004809) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	STOML1 (NM_004809) Human Tagged ORF Clone Lentiviral Particle
Symbol:	STOML1
Synonyms:	hUNC-24; SLP-1; STORP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004809
ORF Size:	1194 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206543).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004809.3 , NP_004800.2
RefSeq Size:	2041 bp
RefSeq ORF:	1197 bp
Locus ID:	9399
UniProt ID:	Q9UBI4
Cytogenetics:	15q24.1
Domains:	Band_7, SCP2
Protein Families:	Transmembrane



[View online »](#)

MW: 43 kDa

Gene Summary: May play a role in cholesterol transfer to late endosomes (PubMed:19696025). May play a role in modulating membrane acid-sensing ion channels. Can specifically inhibit proton-gated current of ASIC1 isoform 1. Can increase inactivation speed of ASIC3. May be involved in regulation of proton sensing in dorsal root ganglions (By similarity). May play a role in protecting FBXW7 isoform 3 from degradation (PubMed:23082202).[UniProtKB/Swiss-Prot Function]