

Product datasheet for **RC206545L1V**

HABP2 (NM_004132) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HABP2 (NM_004132) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HABP2
Synonyms:	FSAP; HABP; HGFAL; NMTC5; PHBP
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_004132
ORF Size:	1680 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206545).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004132.2
RefSeq Size:	3019 bp
RefSeq ORF:	1683 bp
Locus ID:	3026
UniProt ID:	Q14520
Cytogenetics:	10q25.3
Domains:	KR, Tryp_SPc, EGF, EGF
Protein Families:	Druggable Genome, Protease, Secreted Protein



[View online »](#)

MW: 62.7 kDa

Gene Summary: This gene encodes a member of the peptidase S1 family of serine proteases. The encoded preproprotein is secreted by hepatocytes and proteolytically processed to generate heavy and light chains that form the mature heterodimer. Further autoproteolysis leads to smaller, inactive peptides. This extracellular protease binds hyaluronic acid and may play a role in the coagulation and fibrinolysis systems. Mutations in this gene are associated with nonmedullary thyroid cancer and susceptibility to venous thromboembolism. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed. [provided by RefSeq, Jan 2016]