

Product datasheet for RC206557L1

MITF (NM_198158) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: MITF (NM_198158) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: MITF

Synonyms: bHLHe32; CMM8; COMMAD; MI; WS2; WS2A

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC206557).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_198158

ORF Size: 1239 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



MITF (NM_198158) Human Tagged Lenti ORF Clone - RC206557L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 198158.1</u>

RefSeq Size:4454 bpRefSeq ORF:1242 bpLocus ID:4286

 UniProt ID:
 O75030

 Cytogenetics:
 3p13

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Melanogenesis, Melanoma, Pathways in cancer

MW: 46.3 kDa

Gene Summary: The protein encoded by this gene is a transcription factor that contains both basic helix-loop-


helix and leucine zipper structural features. The encoded protein regulates melanocyte development and is responsible for pigment cell-specific transcription of the melanogenesis

enzyme genes. Heterozygous mutations in the this gene cause auditory-pigmentary syndromes, such as Waardenburg syndrome type 2 and Tietz syndrome. [provided by


RefSeq, Aug 2017]



Product images:



Circular map for RC206557L1



Double digestion of RC206557L1 using Sgfl and Mlul