

Product datasheet for **RC206574L2V**

SPINK1 (NM_003122) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SPINK1 (NM_003122) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SPINK1 |
| Synonyms: | PCTT; PSTI; Spink3; TATI; TCP |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_003122 |
| ORF Size: | 237 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC206574). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_003122.2 |
| RefSeq Size: | 438 bp |
| RefSeq ORF: | 240 bp |
| Locus ID: | 6690 |
| UniProt ID: | P00995 |
| Cytogenetics: | 5q32 |
| Protein Families: | Druggable Genome, Secreted Protein |
| MW: | 8.51 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene is a trypsin inhibitor, which is secreted from pancreatic acinar cells into pancreatic juice. It is thought to function in the prevention of trypsin-catalyzed premature activation of zymogens within the pancreas and the pancreatic duct. Mutations in this gene are associated with hereditary pancreatitis and tropical calcific pancreatitis. [provided by RefSeq, Oct 2008]