

Product datasheet for **RC206592L1V**

IDO1 (NM_002164) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IDO1 (NM_002164) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IDO1
Synonyms:	IDO; IDO-1; INDO
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_002164
ORF Size:	1209 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206592).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002164.3
RefSeq Size:	1944 bp
RefSeq ORF:	1212 bp
Locus ID:	3620
UniProt ID:	P14902
Cytogenetics:	8p11.21
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Tryptophan metabolism




[View online »](#)

MW: 45.3 kDa

Gene Summary: This gene encodes indoleamine 2,3-dioxygenase (IDO) - a heme enzyme that catalyzes the first and rate-limiting step in tryptophan catabolism to N-formyl-kynurenine. This enzyme acts on multiple tryptophan substrates including D-tryptophan, L-tryptophan, 5-hydroxy-tryptophan, tryptamine, and serotonin. This enzyme is thought to play a role in a variety of pathophysiological processes such as antimicrobial and antitumor defense, neuropathology, immunoregulation, and antioxidant activity. Through its expression in dendritic cells, monocytes, and macrophages this enzyme modulates T-cell behavior by its peri-cellular catabolization of the essential amino acid tryptophan.[provided by RefSeq, Feb 2011]

Product images:



[RC206592L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC206592L1V particle to overexpress human IDO1-Myc-DDK fusion protein.