

Product datasheet for RC206593L3

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

HSD3B2 (NM_000198) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: HSD3B2 (NM_000198) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: HSD3B2

Synonyms: HSD3B; HSDB; SDR11E2

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC206593).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_000198

ORF Size: 1116 bp





HSD3B2 (NM_000198) Human Tagged Lenti ORF Clone - RC206593L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000198.2</u>

RefSeq Size: 1730 bp RefSeq ORF: 1119 bp

Locus ID: 3284

UniProt ID: P26439

Cytogenetics: 1p12

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Androgen and estrogen metabolism, C21-Steroid hormone metabolism, Metabolic pathways

MW: 42.1 kDa

Gene Summary: The protein encoded by this gene is a bifunctional enzyme that catalyzes the oxidative

conversion of delta(5)-ene-3-beta-hydroxy steroid, and the oxidative conversion of

ketosteroids. It plays a crucial role in the biosynthesis of all classes of hormonal steroids. This gene is predominantly expressed in the adrenals and the gonads. Mutations in this gene are associated with 3-beta-hydroxysteroid dehydrogenase, type II, deficiency. Alternatively spliced

transcript variants have been found for this gene. [provided by RefSeq, Oct 2009]