

## Product datasheet for **RC206613L4V**

### CD1 (CD1D) (NM\_001766) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CD1 (CD1D) (NM_001766) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CD1
Synonyms:	CD1A; R3; R3G1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001766
ORF Size:	1005 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206613).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001766.3</a> , <a href="#">NP_001757.1</a>
RefSeq Size:	3795 bp
RefSeq ORF:	1008 bp
Locus ID:	912
UniProt ID:	<a href="#">P15813</a>
Cytogenetics:	1q23.1
Domains:	IGc1
Protein Families:	Druggable Genome, Transmembrane


[View online »](#)

**Protein Pathways:** Hematopoietic cell lineage

**MW:** 37.7 kDa

**Gene Summary:** This gene encodes a divergent member of the CD1 family of transmembrane glycoproteins, which are structurally related to the major histocompatibility complex (MHC) proteins and form heterodimers with beta-2-microglobulin. The CD1 proteins mediate the presentation of primarily lipid and glycolipid antigens of self or microbial origin to T cells. The human genome contains five CD1 family genes organized in a cluster on chromosome 1. The CD1 family members are thought to differ in their cellular localization and specificity for particular lipid ligands. The protein encoded by this gene localizes to late endosomes and lysosomes via a tyrosine-based motif in the cytoplasmic tail. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016]