

Product datasheet for **RC206629A1V**

Human Interleukin 34 (IL34) (NM_152456) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Interleukin 34 (IL34) (NM_152456) AAV Particle
Tag:	Myc-DDK
Symbol:	Interleukin 34
Synonyms:	C16orf77; IL-34
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206629 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCCGGGGCTTCACCTGGCTGCGCTATCTTGGGATCTTCCTTGGCGTGGCCTTGGGGAATGAGCCTT
TGGAGATGTGGCCCTTGACGCAGAATGAGGAGTGCACTGTACGGGTTTTCTGCGGGACAAGCTGCAGTA
CAGGAGCCGACTTCAGTACATGAACTACTTCCCCATCACTACAAGATCAGTGTGCCTTACGAGGGG
GTGTTTCAGAATCGCCAACGTACCAGGCTGCAGAGGGCCAGGTGAGCGAGCGGGAGCTGCGGTATCTGT
GGTCTTGGTGAGCCTCAGTGCCACTGAGTCGGTGCAGGACGTGCTGCTCGAGGGCCACCATCCTGGAA
GTACCTGCAGGAGGTGCAGACGCTGCTGCTGAATGTCCAGCAGGGCCTCACGGATGTGGAGGTCAAGCCCC
AAGGTGGAATCCGTGTTGTCCCTCTTGAATGCCCGAGGGCCAAACCTGAAGCTGGTGCAGGCCAAAGCCC
TGCTGGACAACTGCTTCCGGGTCATGGAGCTGCTGTACTGCTCCTGCTGTAAACAAAGCTCCGTCTAAA
CTGGCAGGACTGTGAGGTGCCAAGTCTCAGTCTTGACAGCCAGAGCCCTCATTGCAGTATGCGGCCACC
CAGCTGTACCCTCCGCCCCCGTGGTCCCCAGCTCCCCGCCTCACTCCACGGGCTCGGTGAGGCCGTCA
GGGCACAGGGCGAGGGCCTCTTGCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206629 protein sequence Red=Cloning site Green=Tags(s) MPRGFTWLRYLGI FLGVALGNEPLEMWPLTQNEECTVTGFLRDKLQYRSRLQYMKHYFPIN YKISVPYEG VFRIANVTRLQRAQVSERELRYLWVLVSL SATESVQDVLLLEGHP SWKYLQEVQTLLLNVQQGLTDVEVSP KVESVLSLLNAPGPNLKLVRPKALLDNCFRVMELLYCSCCKQSSVLNWQDCEVPSPQSCSPEPSLQYAAT QLYPPPPWSPSPPHSTGSRPVRAQGEGLLP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152456
ORF Size:	726 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152456.1, NP_689669.1</u>
RefSeq Size:	1796 bp
RefSeq ORF:	729 bp
Locus ID:	146433
UniProt ID:	<u>Q6ZMJ4</u>
Cytogenetics:	16q22.1
MW:	27.5 kDa