

Product datasheet for **RC206645L4V**

RHPN1 (NM_052924) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RHPN1 (NM_052924) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RHPN1
Synonyms:	ODF5; RHOPHILIN; RHPN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_052924
ORF Size:	2010 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206645).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_052924.2
RefSeq Size:	3795 bp
RefSeq ORF:	2013 bp
Locus ID:	114822
UniProt ID:	Q8TCX5
Cytogenetics:	8q24.3
Domains:	HR1, PDZ, BRO1
MW:	73.6 kDa


[View online »](#)

Gene Summary:

Has no enzymatic activity. May serve as a target for Rho, and interact with some cytoskeletal component upon Rho binding or relay a Rho signal to other molecules.[UniProtKB/Swiss-Prot Function]