

Product datasheet for **RC206668A1V**

Human RTL6 (NM_032287) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human RTL6 (NM_032287) AAV Particle
Tag:	Myc-DDK
Symbol:	RTL6
Synonyms:	dj1033E15.2; LDOC1L; Mar6; Mart6; SIRH3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206668 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTGCAGCCGAGACGTCCAAAGCTGAAAGCCAGCCTTGGCAGCGTCTCCGAATGCCAGATGGATG
ACGTTATTGACACCCTGACCTCCCTGCGCCTACCAACTCGGCGCTGAGGCGGGAGGCTTCCACCCTGCG
GGCGGAGAAGGCCAATCTCACCACATGCTGGAGAGCGTGATGGCAGAGCTGACCTTGTTACGCACCAGG
GCGCGGATCCCGGGGGCTCTGCAGATCACCCCGCCCATCTCCTCAATTACTTCAAACGGGACTCGACCCA
TGACCACACCTCCAACCTCTCTGCCGAGCCCTTTTCCGGGGACCCAGGCCGGTTGGCGGGTTCTCTGAT
GCAGATGGACAGATTATGATCTTCCAGGCCTCCCGCTTCCCGGGTGAAGCCGAGCGTGTGGCCTTCCTT
GTGTCTCGACTGACTGGGGAGGCGGAGAAGTGGCTATCCCCACATGCAACCTGACAGCCCTTGCGCA
ACAACTATCAGGGGTTCTGGCAGAGTTGCGGAGAACCTACAAGTCTCCGCTCCGGCATGCGCGGCGCGC
CCAAATCAGGAAGACTTCTGCCTCTAATAGGGCTGTGCGAGAGAGGCAGATGCTCTGCCGCCAGCTGGCC
TCTGCGGGCACGGGGCTTGCCAGTGCATCCAGCTTCCAACGGGACTAGTCCAGCGCCAGCCCTGCCTG
CCCGAGCACGGAATCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206668 protein sequence Red=Cloning site Green=Tags(s) MVQPQTSKAESPALAASPNAQMDDVIDTLTSLRLTNSALRREASTLRAEKANLTNMLESVMAELTLLRTR ARIPGALQITPPISSITSNTRPMTTPPTSLPEPFSGDPGRLAGFLMQMDRFMIFQASRFPGEAERVAFL VSRLTGEAEKWAIPHMQPDSPLRNNYQGFLAELRRTYKSPLRHARRAQIRKTSASNRAVRERQMLCRQLA SAGTGPCPVHPASNGTSPAPALPARARNL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_032287
ORF Size:	717 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032287.2</u>
RefSeq Size:	5355 bp
RefSeq ORF:	720 bp
Locus ID:	84247
UniProt ID:	<u>Q6ICC9</u>
Cytogenetics:	22q13.31
MW:	26.2 kDa