

Product datasheet for RC206670

C20orf85 (NM_178456) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C20orf85 (NM_178456) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C20orf85
Synonyms: bA196N14.1; LLC1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC206670 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGCAGAAACCTCTCAGCACCGCGCGGCTGAACGCATGAACCTTGTGGTCAGGATGAGATCTGGA
AATACCGTCTGAAGGCTGAATCGGAAGCACGGCAGAACTGGCCCCAGAACTGGGGTTTTTAACAACCCC
TTTTGAGGAGTTGATCAAGTGTGAAGAAGATCTCCCCACCCAAAGCCAAAATCGAGCTTCTGAGCGT
TTCCGCATCCGGCCGGTACCCAGTGGAGAAGTACATCAAGGTCTTTCCATCCCCCAGTCCCACAGA
CCACCCAGGGCTTCATCGGCTGGAGATCTGCAGTGCCAGGCCTGAACAAGTGCCTAGAGCTCGACGATGC
GATCAGAAGCTGCAAAGGTGCTTTTGCACGAGAGCTGTGCTGGCCAAGCAGGGCGTGCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206670 protein sequence
 Red=Cloning site Green=Tags(s)
 MAQKPLSTAAAERMNLVGQDEIWKYRLKAESARQNWPNWGFLLTPFEELIKCEEDLPTPKPKIELPER
 FRIRPVPVEKYIKVFPSPVPQTQGFIGWRSVPLNKCLELDDAIRSCKGAFARELCWPKQGVH


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6359_g08.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_178456

ORF Size: 411 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178456.3](#)

RefSeq Size: 792 bp

RefSeq ORF: 414 bp


Locus ID: 128602

UniProt ID: [Q9H1P6](#)


Cytogenetics: 20q13.32

MW: 15.7 kDa


Product images:



Circular map for RC206670



Western blot validation of overexpression lysate (Cat# [LY405929]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC206670 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified C20orf85 protein (Cat# [TP306670]). The protein was produced from HEK293T cells transfected with C20orf85 cDNA clone (Cat# RC206670) using MegaTran 2.0 (Cat# [TT210002]).