

Product datasheet for **RC206672A1V**

Human C11orf42 (NM_173525) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C11orf42 (NM_173525) AAV Particle
Tag: Myc-DDK
Symbol: C11orf42
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC206672 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTTGGTGGGTACCCCAACCTGCTGACACTGGATGAAGCTGATGCCACCTGGACCCTCATCAAGGATA
AGGTCATCGAGGAGCACTTTGGGCCAATGCAGTAGCAGTACCTTTCCTGTCAGATGCAGCCTGCTATGA
CCTACTGGGTGTGCTGGTAAACAGTCCCGCCAGCCCATACCCGCTGGCTTTGCCAGGTCGGCAGGGC
CGGAGGGCACTGAAACCAGTGGGGCCACTACCAAGCCTCCTGGAGCAGGCAGGATCTGAGGGTGCCTTCG
CCCACTGCACTCGGGAATACTACCAAATGGCCGAGCAGAGAGAGCCTATGAAGAGACGCGAATGTTGGA
TGGACAGCCCTGCAAGATCCGCCTACATATGGGTGACCTGCGCAAGAAGTTGCCTTCTGTTGCTGCCA
CCAGGGCAGGTGAGCCTACAGCAGACTTTCCTGGCTCCGAAGCACCCACAGCATCTATGTCATCTACC
AGGTCTTCTCTTGTCTGGCTGCAGCTGGGGCTGACGTCTACAGCCCGTGAGCCCCAGCTCCTCCGGCT
ACTTCGGTCATTACCTGTTGCCTTCTCCTGCCTCAAGTTTTCACTGCAGTCTAAGGGCGTGCTGGGACCA
CAGAAGCCTCTACTAAAGACCCATTGCCCCATGGGGCCAACTGGGTGAGACCAACCTCAGCATCATGC
CGCCTCTGGCCCCACATCAGCATCTGCTGATACTGAAGCTGCTGATGTGCCCCACCTGTCCCAGC
CCCACCTACGCCACCTCCCAGGAAGGGCCAGAGGACAAACCCACCAGATTCTCTACAAGGGCCGAAAC
CCCTTCTGGAGGGGGCCAGATACTGTCAGAGAACTGGCTCTTCAGCCCCGAGCCCTCCACCAGGAG
CCCAGGGTGGGGGCCCCAGGGACCCGACGGGCACTCCATGTCCTGCCCTGCTGCAGGGTCTATCCTC
AGAGTTCGACAGTGACGAC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206672 protein sequence Red=Cloning site Green=Tags(s) MLVGTPNLLTLDEADATWTLIKDKVIEEHFGPNAVAVPFLSDAACDYLGLVKQSRPAHTRLALPGRQG RRALKPVGPLPSLLEQAGSEGAFAHCTREYSPNGRAERAYEETRMLDGGQCKIRLHMGDLRKKVAFLLLP PGQVSLQQTLPWLRSTHSIYVIYQVFSWSLQLGLTSTAREPQLRLRLSLPVAFSCLKFSLQSKGVLGP QKPLTKDPLPHGANWVRPNLSIMPPLAPTSASADTTEAADVPPVPAPPTPPPQEGPEDKPTRFSYKGRN PFWRGPIILSENWLFSPRSPPPGAQGGGPRDPDGHMSLPLLQGLSSEFDSD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_173525
ORF Size:	999 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_173525.1, NP_775796.1</u>
RefSeq Size:	1149 bp
RefSeq ORF:	1002 bp
Locus ID:	160298
UniProt ID:	<u>Q8N5U0</u>
Cytogenetics:	11p15.4
MW:	36.4 kDa