

Product datasheet for **RC206673L4V**

STK22C (TSSK3) (NM_052841) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	STK22C (TSSK3) (NM_052841) Human Tagged ORF Clone Lentiviral Particle
Symbol:	STK22C
Synonyms:	SPOGA3; STK22C; STK22D; TSK3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_052841
ORF Size:	804 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206673).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_052841.3
RefSeq Size:	1330 bp
RefSeq ORF:	807 bp
Locus ID:	81629
UniProt ID:	Q96PN8
Cytogenetics:	1p35.1
Protein Families:	Druggable Genome, Protein Kinase
MW:	29.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a kinase expressed exclusively in the testis that is thought to play a role in either germ cell differentiation or mature sperm function. [provided by RefSeq, Jul 2008]