

Product datasheet for **RC206685A1V**

Human Granzyme H (GZMH) (NM_033423) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Granzyme H (GZMH) (NM_033423) AAV Particle
Tag:	Myc-DDK
Symbol:	Granzyme H
Synonyms:	CCP-X; CGL-2; CSP-C; CTLA1; CTSLGL2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206685 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAGCCATTCTCTCTCTGTTGGCCTTTCTTCTGACCCCTGGGGCTGGGACAGAGGAGATCATCGGGG
GCCATGAGGCCAAGCCCCACTCCCGCCCTACATGGCCTTTGTTCAAGTTCTGCAAGAGAAGAGTCGGAA
GAGGTGTGGCGGCATCCTAGTGAGAAAGGACTTTGTGCTGACAGCTGCTCACTGCCAGGGAAGCTCCATA
AATGTCACCTTTGGGGGCCACAATATCAAGGAACAGGAGCGGACCCAGCAGTTTATCCCTGTGAAAAGAC
CCATCCCCCATCCAGCCTATAATCCTAAGAACTTCTCCAACGACATCATGCTACTGCAGCTGGAGAGAAA
GGCCAAGTGGACCACAGCTGTGCGGCCCTCTCAGGCTACCTAGCAGCAAGGCCAGGTGAAGCCAGGGCAG
CTGTGCACTGTGGCTGGCTGGGTTATGTCTCAATGAGCACTTTAGCAACCACACTGCAGGAAGTGTTC
TGACAGTGCAGAAGGACTGCCAGTGTGAACGTCTCTTCATGGCAATTACAGCAGAGCCACTGAGATTTG
TGTGGGGGATCCAAAGAAGACACAGACCGGTTTCAAGGGGACTCCGGGGGGCCCTCGTGTGTAAGGAC
GTAGCCCAAGGTATTCTCTCTATGAAACAAAAAGGGACACCTCCAGGAGTCTACATCAAGGTCTCAC
ACTTCTGCCCTGGATAAAGAGAACAATGAAGCGCCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206685 protein sequence Red=Cloning site Green=Tags(s) MQPFLLLLAFLLTPGAGTEEIIGGHEAKPHSRPYMAFVQFLQEKSRRKCGGILVRKDFVLTAAHCQGSSI NVTLGAHNIKEQERTQQFIPVKRPIPHPAYNPKNFSNDIMLLQLERKAKWTTAVRPLRLPSSKAQVKPGQ LCSVAGWGYVSMSTLATTLEVLLTVQKDCQERLFHGNYSRATEICVGDPKKTQTGFKGDSGGPLVCKD VAQGILSYGNKKGTPPGVYIKVSHFLPWIKRTMKRL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_033423
ORF Size:	738 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_033423.2</u>
RefSeq Size:	987 bp
RefSeq ORF:	741 bp
Locus ID:	2999
UniProt ID:	<u>P20718</u>
Cytogenetics:	14q12
MW:	27.3 kDa