

Product datasheet for **RC206691L2V**

KLF4 (NM_004235) Human Tagged ORF Clone Lentiviral Particle


Product data:

Product Type:	Lentiviral Particles
Product Name:	KLF4 (NM_004235) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KLF4
Synonyms:	EZF; GKLF
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004235
ORF Size:	1410 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206691).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004235.3 , NP_004226.2
RefSeq Size:	2639 bp
RefSeq ORF:	1440 bp
Locus ID:	9314
UniProt ID:	O43474
Cytogenetics:	9q31.2
Domains:	zf-C2H2



[View online »](#)

Protein Families:	Adult stem cells, Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Transcription Factors
MW:	49.9 kDa
Gene Summary:	This gene encodes a protein that belongs to the Kruppel family of transcription factors. The encoded zinc finger protein is required for normal development of the barrier function of skin. The encoded protein is thought to control the G1-to-S transition of the cell cycle following DNA damage by mediating the tumor suppressor gene p53. Mice lacking this gene have a normal appearance but lose weight rapidly, and die shortly after birth due to fluid evaporation resulting from compromised epidermal barrier function. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]

Product images:

[RC206691L2] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC206691L2V particle to overexpress human KLF4-mGFP fusion protein.