

#### Product datasheet for RC206692L1

### https://www.origene.com techsupport@origene.com EU: info-de@origene.com

# CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

## PTGR1 (NM\_012212) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type: Expression Plasmids** 

**Product Name:** PTGR1 (NM\_012212) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: PTGR1

Synonyms: DIG-1; LTB4DH; PGR1; ZADH3

**Mammalian Cell** 

Selection:

None

Vector: pLenti-C-Myc-DDK (PS100064) E. coli Selection: Chloramphenicol (34 ug/mL)


**ORF Nucleotide** 

The ORF insert of this clone is exactly the same as(RC206692).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

ACCN: NM\_012212

**ORF Size:** 987 bp



#### PTGR1 (NM\_012212) Human Tagged Lenti ORF Clone - RC206692L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 012212.2</u>, <u>NP 036344.1</u>

 RefSeq Size:
 1249 bp

 RefSeq ORF:
 990 bp

 Locus ID:
 22949

 UniProt ID:
 Q14914

 Cytogenetics:
 9q31.3

**Protein Families:** Druggable Genome

**MW:** 35.9 kDa

**Gene Summary:** This gene encodes an enzyme that is involved in the inactivation of the chemotactic factor,


leukotriene B4. The encoded protein specifically catalyzes the NADP+ dependent conversion

of leukotriene B4 to 12-oxo-leukotriene B4. A pseudogene of this gene is found on chromosome 1. Alternative splicing results in multiple transcript variants. [provided by

RefSeg, Mar 2009]



# **Product images:**



Circular map for RC206692L1