

## Product datasheet for **RC206697L3V**

### C4orf22 (CFAP299) (NM\_152770) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	C4orf22 (CFAP299) (NM_152770) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CFAP299
Synonyms:	C4orf22
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_152770
ORF Size:	699 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206697).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152770.1</a> , <a href="#">NP_689983.1</a>
RefSeq Size:	903 bp
RefSeq ORF:	702 bp
Locus ID:	255119
UniProt ID:	<a href="#">Q6V702</a>
Cytogenetics:	4q21.21
MW:	26.9 kDa
Gene Summary:	May be involved in spermatogenesis.[UniProtKB/Swiss-Prot Function]


[View online »](#)