

Product datasheet for **RC206701L4V**

ANKRD53 (NM_024933) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ANKRD53 (NM_024933) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ANKRD53
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024933
ORF Size:	1029 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206701).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024933.2 , NP_079209.2
RefSeq Size:	2120 bp
RefSeq ORF:	1032 bp
Locus ID:	79998
UniProt ID:	Q8N9V6
Cytogenetics:	2p13.3
MW:	38.2 kDa
Gene Summary:	Required for normal progression through mitosis. Involved in chromosome alignment and cytokinesis via regulation of microtubules polymerization.[UniProtKB/Swiss-Prot Function]


[View online »](#)