

Product datasheet for **RC206706L4V**

KLHDC8B (NM_173546) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KLHDC8B (NM_173546) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KLHDC8B
Synonyms:	CHL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_173546
ORF Size:	1062 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206706).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_173546.2
RefSeq Size:	1999 bp
RefSeq ORF:	1065 bp
Locus ID:	200942
UniProt ID:	Q8IXV7
Cytogenetics:	3p21.31
MW:	38.1 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein which forms a distinct beta-propeller protein structure of kelch domains allowing for protein-protein interactions. Mutations in this gene have been associated with Hodgkin lymphoma. [provided by RefSeq, Sep 2010]