

Product datasheet for **RC206715A1V**

Human DHH (NM_021044) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DHH (NM_021044) AAV Particle
Tag:	Myc-DDK
Symbol:	DHH
Synonyms:	GDMN; GDXYM; HHG-3; SRXY7
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206715 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTCTCTGACCAATCTACTGCCCTGTGCTGCTTGGCACTTCTGGCGTGCCAGCCAGAGCTGCG
GGCCGGGCGGGGCGGTTGGCCGGCGCCGCTATGCGCGCAAGCAGCTCGTGCCGCTACTCTACAAGCA
ATTTGTGCCGGCGTGCCAGAGCGGACCTGGGCGCCAGTGGGCCAGCGAGGGGAGGGTGGCAAGGGGC
TCCGAGCGCTTCCGGGACCTCGTGCCCAACTACAACCCGACATCATCTTCAAGGATGAGGAGAACAGTG
GAGCCGACCGCTGATGACCGAGCGTTGTAAGGAGCGGGTGAACGCTTTGGCCATTGCCGTGATGAACAT
GTGGCCCGGAGTGCGCTACGAGTGACTGAGGGCTGGGACGAGGACGGCCACCACGCTCAGGATTCACTC
CACTACGAAGGCCGTGCTTTGGACATCACTACGTCTGACCGCGACCGCAACAAGTATGGGTTGCTGGCGC
GCCTCGCAGTGGAAGCCGGCTTCGACTGGGTCTACTACGAGTCCCGCAACCACGTCCACGTGTCGGTCAA
AGCTGATAACTCACTGGCGGTCCGGGCGGGCGGCTGCTTTCCGGGAAATGCAACTGTGCGCCTGTGGAGC
GGCGAGCGGAAAGGGCTGCGGGAATGCACCGCGGAGACTGGTTTTGGCGGCCGATGCGTCAGGCCGGG
TGGTGCCACGCCGTGCTGCTTCTTGACCGGGACTTGACGCGCCGGGCTTCATTTGGCTGTGGA
GACCGAGTGGCTCCACGAACTGTTGCTCACGCCCTGGCACCTGGTGTTCGCGCTCAGGGCCGGCG
CCGCGCCAGGCGACTTTGCACCGGTGTTGCGCGCCGGCTACGCGCTGGGACTCGGTGCTGGCGCCCG
GCGGGGATGCGCTTCGGCCAGCGCGCTGGCCCGTGTGGCGCGGGAGGAAGCCGTGGGCGTGTTGCGGCC
GCTCACCGCGCACGGGACGCTGCTGGTGAACGATGTCCTGGCCTCTTGCTACGCGGTTCTGGAGAGTCAC
CAGTGGGCGCACCGCGCTTTTGGCCCTTGAGACTGCTGCACGCGCTAGGGGCGCTGCTCCCGCGGGG
CCGTCCAGCCGACTGGCATGCATTGGTACTCTCGGCTCTCTACCGCTTAGCGGAGGAGCTACTGGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206715 protein sequence Red=Cloning site Green=Tags(s) MALLTNLLPLCCLALLALPAQSCGPGRGPVGRRRYARKQLVPLLYKQFVPGVPERTLGASGPAEGRVARG SERFRDLVPNYPNDIIFKDEENSGADRLMTERCKERVNALAIIVMMNWPVRLRVTEGWDEGDGHAQDSL HYEGRALDITTSDDRNRKYGLLARLAVEAGFDWVYYESRNHVVHVSVKADNSLAVRAGGCFPGNATVRLWS GERKGLRELHRGDWVLAADASGRVPTPVLLFLDRDLQRRASFVAVETEWPPRLLLLTPWHLVFAARGPA PAPGDFAPVFAARRLRAGDSVLAPGGDALRPARVARVAREEAVGVFAPLTAHGTLVNDVLASCYAVLESH QWAHRAFAPLRLLHALGALLPGGAVQPTGMHWYSRLLYRLAEELLG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_021044
ORF Size:	1188 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_021044.2
RefSeq Size:	1971 bp
RefSeq ORF:	1191 bp
Locus ID:	50846
UniProt ID:	O43323
Cytogenetics:	12q13.12
MW:	43.6 kDa