

Product datasheet for **RC206725L3V**

DMRTC2 (NM_001040283) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DMRTC2 (NM_001040283) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DMRTC2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001040283
ORF Size:	1101 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206725).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001040283.1
RefSeq Size:	1821 bp
RefSeq ORF:	1104 bp
Locus ID:	63946
UniProt ID:	Q8IXT2
Cytogenetics:	19q13.2
Protein Families:	Transcription Factors
MW:	39.1 kDa
Gene Summary:	May be involved in sexual development.[UniProtKB/Swiss-Prot Function]



[View online »](#)