

Product datasheet for **RC206735L4V**

IGSF11 (NM_152538) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IGSF11 (NM_152538) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IGSF11
Synonyms:	BT-IgSF; CT119; CXADRL1; Igsf13; VSIG3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152538
ORF Size:	1290 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206735).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152538.2
RefSeq Size:	3577 bp
RefSeq ORF:	1293 bp
Locus ID:	152404
UniProt ID:	Q5DX21
Cytogenetics:	3q13.32
Domains:	ig, IGv, IGc2, IG
Protein Families:	Transmembrane



[View online »](#)

MW: 46.2 kDa

Gene Summary: IGSF11 is an immunoglobulin (Ig) superfamily member that is preferentially expressed in brain and testis. It shares significant homology with coxsackievirus and adenovirus receptor (CXADR; MIM 602621) and endothelial cell-selective adhesion molecule (ESAM).[supplied by OMIM, Apr 2005]