

Product datasheet for **RC206736L3V**

PTPNS1L2 (SIRPD) (NM_178460) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PTPNS1L2 (SIRPD) (NM_178460) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PTPNS1L2
Synonyms:	dj576H24.4; PTPNS1L2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178460
ORF Size:	591 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206736).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178460.1 , NP_848555.1
RefSeq Size:	868 bp
RefSeq ORF:	594 bp
Locus ID:	128646
UniProt ID:	Q9H106
Cytogenetics:	20p13
Protein Families:	Druggable Genome, Phosphatase
MW:	21.7 kDa



[View online »](#)